NTI-Kahla GmbH • Rotary Dental Instruments
Im Camisch 3 • D-07768 Kahla/Germany
Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
e-mail: nti@nti.de • http://www.nti.de

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
GB Product Information Sheet - Tissue Trimmer

Description of the product:

The Tissue Trimmer is a surgical instrument for use on the gingiva. The tip of this instrument is made of special ceramic.

Range of indications:

- Widening the sulcus for impression taking and periodontal treatment
- Exposing deep cervical cavities
- Recovering impacted teeth and intraosseous implants
- Papillectomies
- Contouring the gingiva

Application:

The Tissue Trimmer must only be used in a turbine. The working speed is 300,000 - 500,000 r.p.m. Lower speeds may cause injury and must be avoided.

Coolant spray:

The turbine coolant spray must be switched off. Coolant spray impedes coagulation in the capillaries which may lead to increased bleeding.

Disinfection:

The ceramic tip is more delicate than tungsten carbide or diamond-coated instruments. Therefore, please do not clean this instrument in an ultrasonic cleaner with other instruments.

Please use only approved disinfectants suitable for use with rotary instruments. The disinfectant must comply with legal requirements of and be licensed in the land in question. Alcohol is not suitable for use and does not provide sufficient disinfection.

Please comply with the manufacturer's regulations for procedures before and after disinfecting.

Sterilisation:

Use an autoclave for sterilisation.

Important:

Throughout the whole preparation procedure ensure that the tip of the instrument is not damaged. Check for damage after every disinfection process using a magnifying-glass of 10x magnification.

Detailed recommendations for preparation in accordance with DIN EN ISO 17664 www.nti.de

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